







# Instruction manual GMS\_2CH\_HP Gas mixer

Intended use .....	1
Technical Data .....	2
Installation .....	2
Driver installation for USB - RS485 converter - connection to PC .....	2
Software – Installation .....	2
Installation of GMS_2CH ( hardware ) .....	3
1. Electrical connection .....	3
2. Gas connections.....	3
Startup operations.....	3
Decommission.....	3
Maintenance.....	4
Calibration .....	4
EU Declaration of Conformity .....	5

## **Intended use**

<p>Before operation this instruction manual has to be read. During installation, operation, maintenance and shut-off this instruction manual has to be observed.</p> <p>The GMS_2CH is designed for use of the following gases :</p> <ul style="list-style-type: none"> <li>• Air</li> <li>• Nitrogen</li> <li>• Oxygen</li> <li>• Carbon Dioxide</li> <li>• other non-corrosive gases</li> </ul> <p><b>All used gases must be free of particles.</b>  <b>If necessary a particulate filter has to be installed at the gas inlets.</b></p>	
--	---

<p>When the GMS_2CH is operated with toxic, ignitable, combustible or explosive substances all safety instruction for these substances have to be observed carefully.</p> <p>In the event of an accident :</p> <ol style="list-style-type: none"> <li>1. turn OFF the power switch.</li> <li>2. turn OFF all power switches of peripheral units.</li> <li>3. close all main valves of connected gases.</li> <li>4. unplug the power cable.</li> </ol>	   
---	--

<p>When the GMS_2CH is operated with substances, not listed above, the materials compatibility has to be checked for all sample wetted components.</p>	
--	---

QCAL Messtechnik GmbH  
 Gardinistr. 143  
 81375 München  
 GERMANY  
 Email : [info@qcal.de](mailto:info@qcal.de)  
[www.qcal.de](http://www.qcal.de)

### **Important notice :**

For gases the standard volume relates to 0 °C and 1013.25 hPa.  
 Abbreviations and symbols :

MFC = thermal massflow controller for gases



WARNING - important notice



information

### **Technical Data**

Mains power :	88 ~ 264 V AC, 47 ~ 63 Hz
Current consumption :	max. 3.15 A, fused ( 5 x 20 mm )
Ambient conditions :	15 ° C to 40°C, 0 to 95% rH
Gas inlet conditions :	particle-free gases, max. inlet pressure : 8 bar(g)
Dimensions :	250 x 250 x 128 mm
Weight :	3 kg

### **Installation**

#### **Driver installation for USB - RS485 converter - connection to PC**



For WINDOWS 64bit operating systems the driver is installed automatically from the setup program.

In all other cases : before connection of GMS\_2CH with the PC read the instructions manual from the directory /driver. Connect the USB cable and install the driver according to the WINDOWS operating system.




Fig.1: USB connection GMS\_2CH - PC

After connection GMS\_2CH to PC a new COM ( serial ) port will be installed. This COM - port number - shown in the WINDOWS device manager - has to be registered in the GMS\_2CH software.

#### **Software – Installation**

The GMS\_2CH software has to be installed by executing setup.exe from the /setup directory.



For operation MS EXCEL® 2016 must be installed on the computer.



If the error message : **RegSvr32 – Error loading the module 'files'** is displayed, it can be ignored.

## Installation of GMS\_2CH ( hardware )

### 1. Electrical connection

Connect the mains plug to mains voltage : 88 ~ 264 V AC, 47 ~ 63 Hz.

The GMS\_2CH is switched ON / OFF by the red mains switch.



Fig. 2: power switch

Fuse : 3.15 A  
For exchange : **Unplug the mains plug !**

A green LED indicates the correct power supply voltage ( 24V DC ).

### 2. Gas connections



For SWAGELOK fittings : For assembly of tube fittings the Swagelok instructions manual has to be observed.



For installation of tubings the bulkhead union has to be fixed according to Fig.3 with a suitable wrench.



Fig. 3: fixing



Maximum inlet pressure : 10 bar(g)



Maximum outlet pressure : 8 bar(g).



**When the operation is finished : close all valves at the gas containers !**

### Startup operations

1. Switch GMS\_2CH ON - red mains switch. A LED indicates correct power connection.
2. Connect the GMS\_2CH to the PC : USB 2.0 or 3.0 port.
3. Start the GMS\_2CH software.

### Decommission

If the GMS\_2CH is not operated for a long time the valves on the gas cylinders must be closed.

## Maintenance

The enclosure can be cleaned with a mild detergent.



Avoid any moisture or particles inside the gas tubings - that can lead to serious damage of the MFCs.

## Calibration

The (re)calibration can be performed by the user. A reliable reference device is necessary for every gas type.



The calibration data are stored in the MS EXCEL file `GMS_2CH.xlsx`. For each MFC and for each gas exists a worksheet, that can be edited by the user.

Only the yellow fields are editable – the others are locked. They contain the measurement data from each MFC and the corresponding reference data :

Flow\_ref/Ncm<sup>3</sup>/Minute : gas flow of the reference instrument / reference method

Flow\_GMS/Ncm<sup>3</sup>/Minute : gas flow from the `GMS_2CH`.

32 calibration points for each gas are required.



All data from the MFCs refer to 0°C and 1013,25 hPa, so the reference data must also refer to these conditions.

The screenshot shows an Excel spreadsheet with the following data:

	A	B	C	D	E	F	G
1		GAS : air			Flow = slope_i x Flow_GMS + zero_i		
2		ascending values !					
3		Flow_ref/Ncm <sup>3</sup> /Minute	Flow_GMS/Ncm <sup>3</sup> /Minute	slope	zero		reference
4	1	0,000	0,000	1	0	1	Zirox
5	2	0,200	0,200	1	0	2	Zirox
6	3	0,400	0,400	1	0	3	Zirox
7	4	0,600	0,600	1	0	4	Zirox

Fig.4 : Calibration data `GMS_2CH`

When the `GMS_2CH` is started again the program will automatically read the new calibration data.

Issue : February 2024

QCAL Messtechnik GmbH  
Guardinstr. 143  
81375 München  
Germany

Phone : ++49 (0)89 84060347  
email : info@qcal.de

# EU Declaration of Conformity



We



QCAL Messtechnik GmbH  
Guardinstr. 143  
D-81375 München  
GERMANY  
Phone : ++49 (0) 89 84060347  
Fax : ++49 (0) 89 999 64416  
eMail : [info@qcal.de](mailto:info@qcal.de)  
<http://www.qcal.de/>

declare under our sole responsibility that the following products :

Equipment: Gas mixing units  
Model/type: GMS\_2CH\_HP , GMS\_3CH\_HP, GMS\_4CH\_HP

are in conformity with the following harmonised standards

Machinery Directive 2006/42 / EG

Electromagnetic compatibility (EMC) legislation , Part 4-2, 4-3, 4-4, 4-5

Munich, Oct / 09 / 2020

QCAL Messtechnik GmbH

A handwritten signature in blue ink, appearing to read 'Frank Monshausen', is written over a horizontal line.

Dr. Frank Monshausen, Managing Director